

## UbiqVue Mobile Device and OS Compatibility Guide

Updated June 01, 2026

Popular mobile devices and operating systems (OS) are routinely tested to ensure reliable connectivity and compatibility with UbiqVue Biosensors. LifeSignals recommends reviewing this guide before updating your phone's OS or using an app on a new device

This guide contains the most up-to-date compatibility information and is regularly updated as new devices and OS versions are evaluated. Beta OS versions are not supported.

### UbiqVue Holter App Compatibility

**Model:** UA2251-R

**For use with:** UbiqVue 2AYe Holter System

**App versions:** Android 1000002114H | iOS 1000002115J

Device	OS
<ul style="list-style-type: none"><li>iPhone: 11, 13, 13 Mini, SE2, 14 Plus, 12 Mini, 15, 16, 17</li></ul>	<b>iOS:</b> 26.5
<ul style="list-style-type: none"><li>Samsung Galaxy: A06, A15, M16, A25</li><li>Samsung Galaxy: A07, M36 5G</li><li>Xiaomi: Mi 11X</li><li>VIVO: V21 5G</li><li>Pixel 6a</li><li>Moto G67 Power 5G</li><li>Nothing Phone (3a) Lite</li><li>OnePlus Nord CE5</li></ul>	<b>Android:</b> 13, 14, 15, 16