

UbiqVue™ 2AYe Holter System

Featuring 2-Channel ECG Wearable Biosensor with Wi-Fi & Bluetooth® Connectivity



A single-use, chest-worn Holter device that enables 2-channel ECG and event data to be continuously collected and transferred to a secure server via the Patient Holter App or Relay Device ready for storage, analysis and report generation.

Disposable & Affordable

Water resistant, single use Biosensors with up to eight day¹ wear duration

Human Factor Design

Minimal skin preparation required, high patient comfort, convenience and mobility

Reliable and Secure

HIPAA & GDPR compliant encrypted data transfer

Clinical-grade Accuracy

Proven reliability, accuracy, and usability that meet international regulatory requirements

Key features:



2-channel ECG¹



Beat & Rhythm Analysis¹



Heart Rate¹



Heart Rate Variability¹



Event Recorder



Up to 8 day wear duration

UbiqVue™ 2AYe Holter System Product Key Specifications

Biosensor

Product Codes	System - UX2251. Biosensor - UB2251, UB2251A.	
Wireless & Security	WLAN (IEEE 802.11b) 2.400-2.4835 GHz WPA2 - PSK/ CCMP security 5m (typical) CCK and DSSS	BLE 5.2 2.400-2.4835 GHz AES-CCM 128-bit Encryption 5m (typical) FSK
Physical	116 mm x 91 mm x 15.5 mm 32g IP24 water ingress protection 8 days memory cache Patient activated event recorder	
Battery	Lithium Manganese dioxide Li-MnO ₂ Up to 8 days battery life Defibrillation-proof type CF applied part	
Usage	Single use, disposable for home and healthcare settings	

Cloud-Processed Clinical Data

ECG	2 channels 16 bit resolution Monitoring mode : 0.2 Hz to 40 Hz (with 244.14 samples per second)
Beat & Rhythm Analysis	Ventricular Beats, Supraventricular Ectopic Beats, Atrial Fibrillation / Flutter Duration (Burden)
Heart Rate	30-300 bpm
Heart Rate Variability	RMSSD, SDNN

For more information visit www.lifesignals.com



LifeSignals, Inc., 426 S Hillview Drive,
Milpitas, CA 95035, USA



Renew Health Ltd, IDA Business & Technology Park,
Dublin Road, Athlone Co. Westmeath, Ireland



UX2251

1. Cloud server processed from Biosensor biodata.
Consult your sales representative for local market clearance and availability.

LifeSignals, Inc. does not make any claims; for informational purpose only.

LifeSignals, Inc. ©

Doc: MKT-SM-033F August 2026